Product Names: **CRAFT™** and **CRAFT™ DUO-VAC™** Suction Pumps

## ZDOCK105\_EN Rev. 4I

Section	Language	
1	en	CRAFT™ and CRAFT™ DUO-VAC™ Suction Pumps Instructions for Use
2	cn	CRAFT™ 和 CRAFT™ DUO-VAC™ 抽吸泵使用说明
3	de	CRAFT™ und CRAFT™ DUO-VAC™ Saugpumpen Gebrauchsanweisung
4	dk	CRAFT™ og CRAFT™ DUO-VAC™ sugepumper Brugsanvisning
5	es	Bombas de succión CRAFT™ y CRAFT™ DUO-VAC™ Instrucciones De Uso
6	fr	Pompes d'aspiration CRAFT™ et CRAFT™ DUO-VAC™ Mode D'Emploi
7	it	Pompe di aspirazione CRAFT™ e CRAFT™ DUO-VAC™ Istruzioni Per L'Uso
8	nl	CRAFT™ en CRAFT™ DUO-VAC™ Afzuigpompn Gebruiksaanwijzing
9	pl	Pompy ssące CRAFT™ i CRAFT™ DUO-VAC™ Instrukcja Obsługi
10	se	CRAFT™ Och CRAFT™ DUO-VAC™ Sugpumpar Bruksanvisning

1639

**ROCKET MEDICAL PLC**  
Sedling Road, Washington,  
NE38 9BZ, UK

**ROCKET MEDICAL GmbH**  
Am Rosengarten 48,  
15566 Schöneiche.  
Germany

**ROCKET MEDICAL GmbH**  
Am Rosengarten 48,  
15566 Schöneiche.  
Germany

**ROCKET MEDICAL LLC**  
50 Corporate Park Drive #890.  
East Pembroke. MA. 02359.  
USA

**Rocket Medical Pty Ltd**  
Suite 6, 157 Gordon Street,  
Port Macquarie.  
NSW 2444. Australia

Australian Sponsor:

### COPYRIGHT

This manual contains information that is subject to copyright. All rights reserved.

This manual should not be photocopied, duplicated or distributed completely or in part without the approval of Rocket Medical plc

**Product Names: CRAFT™ and CRAFT™ DUO-VAC™ Suction Pumps**

**CONTENTS**

<b>1. SAFETY INSTRUCTIONS</b> .....	1
<b>2. ABOUT CRAFT™ and CRAFT™ DUO-VAC™ PUMPS</b> .....	4
2.2 INTENDED USE.....	6
2.3 INDICATIONS.....	6
2.4 CONTRAINDICATIONS.....	6
2.5 WARNINGS:.....	6
2.6 BENEFIT RISK.....	7
2.7 PERFORMANCE CLAIMS AND CLINICAL BENEFIT.....	7
2.8 UNDESIRABLE SIDE EFFECTS.....	7
2.9 DESCRIPTION OF THE DEVICE:.....	7
2.10 USER.....	7
2.11 PATIENT POPULATION.....	7
2.12 ENVIRONMENT.....	7
<b>3. GENERAL INFORMATION</b> .....	8
3.1 SUPPLY VOLTAGE SELECTION.....	8
3.2 ELECTROMAGNETIC COMPATIBILITY.....	8
3.3 ESSENTIAL CONSUMABLES.....	9
3.4 PACKAGING.....	9
3.5 SUCTION RANGE.....	9
3.6 CALIBRATION.....	9
<b>4. SYMBOLS AND SAFETY SIGNS USED ON THE DEVICE, LABELLING AND ACCESSORIES</b> .....	9
<b>5. OPERATING INSTRUCTIONS</b> .....	11
<b>6. TROUBLESHOOTING</b> .....	13
6.1 BLOCKED VACUUM LINES.....	13
6.2 WATER TRAP SET.....	14
<b>7. DEVICE MAINTENANCE</b> .....	14
<b>8. SERVICE INFORMATION</b> .....	15
8.1 RETURNING THE PUMP.....	15
8.2 SERVICE AGENTS.....	15
<b>9. WARRANTY</b> .....	16
<b>10. TECHNICAL SPECIFICATIONS</b> .....	17
10.1 CLASSIFICATION.....	17
10.2 SPECIFICATIONS.....	17
10.3 EMC TABLES:.....	17
<b>11. DISPOSAL</b> .....	21

*In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at [pncf@rocketmedical.com](mailto:pncf@rocketmedical.com) as well as to the competent health authority in the country that the user/patient resides.*

## I. SAFETY INSTRUCTIONS

This manual describes the operation and intended use of the device and the associated consumables and it is essential that you use this document to familiarise yourself with the correct function and operation of the device before use. Failure to follow these instructions may result in serious injury to the patient or operator and can lead to damage or breakdown of the device. In case the device fails during an operation, a replacement device and replacement disposables should be kept within reach so that the operation may be completed.

This manual does not provide a detailed description of the oocyte harvesting procedure and is not intended as a training guide for users inexperienced in the technique.

The device must be used with R57685 Water Trap Set for CRAFT™ Suction Pumps. The usage of non-approved tubing or filter sets may impair pump performance, lead to increased risk to patients and harvested oocytes and will invalidate the Warranty.

The water trap sets are designed to prevent fluid contamination of the vacuum pump. If the device has been used with a non-approved filter set or there is any evidence or suspicion that the pump may have been contaminated with fluid during use, it must be removed from service and returned for examination immediately. Please contact your nearest Service Centre for advice.

The use of high vacuum levels may lead excessive fluid flow rates which may result in damage to the oocyte and reduced fertilisation rates. Damage to oocytes in harvesting systems is principally caused by turbulent flow which can lead to physical shearing stresses on the cumulus sufficient to denude or damage the fragile zona (Reeves et al 1989). Flow rate is a function needle set configuration and vacuum applied.

HIGH vacuum setting must only be used to clear blockages from a needle set.

Refer to Table I for recommended vacuum setting for given needle set configurations. Refer all servicing to the manufacturer's authorised service agent.

Do not use in an area where flammable gases are present.

<p><b>WARNING</b></p> 	<p><b>READ THIS MANUAL CAREFULLY:</b> Please familiarise yourself with the contents of this manual before attempting to use the device. Failure to observe these instructions may result in damage to the pump or cause injury to the patient or user. This device should only be used by suitably qualified personnel.</p>
<p><b>WARNING</b></p> 	<p><b>Electric Shock hazard:</b> The equipment is to be used only with electrical systems complying with all IEC, CEC and NEC requirements</p>
<p><b>WARNING</b></p> 	<p><b>Any adjustment, modification or repairs to the equipment should be carried out by authorised service agents.</b> <b>This device should only be operated by appropriately qualified personnel</b> <b>Ensure the set vacuum level is appropriate.</b> <b>No user serviceable parts inside</b></p>
<p><b>DISPOSAL</b></p> 	<p><b>Disposal of this device must be undertaken with regard to the WEEE directive (2002/96/EC).</b></p>
<p><b>WARNING</b></p> 	<p>This device must only be operated with approved consumable</p>

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at [pncf@rocketmedical.com](mailto:pncf@rocketmedical.com) as well as to the competent health authority in the country that the user/patient resides.

## Suction Pumps

## INSTRUCTIONS FOR USE

<b>WARNING</b> 	Protect the device from ingress of liquid. Should any liquid enter the device, discontinue use immediately and refer to an authorised service agent
<b>WARNING</b> 	Device can cause explosion in the presence of flammable gases

*In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at [pncf@rocketmedical.com](mailto:pncf@rocketmedical.com) as well as to the competent health authority in the country that the user/patient resides.*

## I. ABOUT CRAFT™ and CRAFT™ DUO-VAC™ PUMPS

**Product Names:** CRAFT™ and CRAFT™ DUO-VAC™ Suction Pumps

**Product Codes:** CRAFT™ Suction Pumps R29655 (220-240V AC 50Hz) and R29654 (110-120V AC 60Hz)  
CRAFT™ DUO-VAC™ Suction Pumps R29660 (220-240V AC 50Hz) and R29661 (110-120V AC 60Hz)

**Device Image:**

**Figure I CRAFT™ SUCTION PUMP**

1. Illuminated I/O (On/Off) power switch
2. Vacuum display (in mmHg)
3. High vacuum (-440 mmHg) control button
4. Vacuum control – clockwise to increase, anticlockwise to decrease the set value
5. Water trap connection port for use with R57685 Water Trap Set for CRAFT™ Suction Pumps
6. R57685 Water Trap Set for CRAFT™ Suction Pumps, single-use only
7. Footswitch connection port
8. Footswitch: medium vacuum 0 to -200 mmHg

*In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at [pncf@rocketmedical.com](mailto:pncf@rocketmedical.com) as well as to the competent health authority in the country that the user/patient resides.*

Figure 2 CRAFT™ DUO-VAC™ Suction Pump

1. Illuminated I/O (On/Off) power switch
2. Vacuum control dial – clockwise to increase, anticlockwise to decrease the set value
3. Vacuum display (in mmHg)
4. Footswitch connection ports
5. Water trap connection port for use with R57685 Water Trap Set for CRAFT™ Suction Pumps
6. R57685 Water Trap Set for CRAFT™ Suction Pumps, single-use only
7. Dual footswitch: medium vacuum 0 to -200 mmHg (black) & high vacuum -440 mmHg (white)

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at [pncf@rocketmedical.com](mailto:pncf@rocketmedical.com) as well as to the competent health authority in the country that the user/patient resides.

## Suction Pumps

## INSTRUCTIONS FOR USE

**2.2 INTENDED USE:** For the generation of medium vacuum between 0 to -200 mmHg to permit the aspiration of follicular fluid, oocytes and ovarian fluid.

**2.3 INDICATIONS:** For use as part of the treatment of infertility relating to in-vitro fertilisation (IVF) and other gynaecological procedures. This device may only be used by, or under the supervision of, appropriately trained personnel in conjunction with clinical practice guidelines.

**2.4 CONTRAINDICATIONS:**

The patient should be made aware of;-

- Not intended for use where ovarian aspiration or the aspiration of ovarian fluid is contraindicated:
- Patients with vaginal infections and or sexually transmitted diseases as this increases the risk of transmission of infection from the pump to the next patient
- Suspected pregnancy as this could result in miscarriage or fetal injury
- Where the pelvic anatomy cannot be accurately assessed by ultrasonography as this could result in damage to other organs and vessels
- Not intended for surgical suction or other applications as the diameter of the tubing is too narrow for adequate drainage.

**1.5 WARNINGS:**

**WARNING** Users should be familiar with and adhere to all warnings, cautions and instructions for use that are labelled on the device and included in the User Manual.

**WARNING** This device should only be used by, or under the supervision of, appropriately trained personnel and in conjunction with current local clinical practice guidelines.

**WARNING** To avoid the risk of electric shock this equipment must only be connected to a supply main with a protective earth.

**WARNING** The device can be isolated from the mains supply by removal of the mains supply cord from the rear of the equipment.

**WARNING** The device is not approved for connection to any other device except the R57685 Water Trap Set for CRAFT™ Suction Pumps

**WARNING** The R57685 Water Trap Set for CRAFT™ Suction Pumps is single use and its instructions for use must be followed at all times.

**WARNING** Do NOT use the R57685 Water Trap Set for CRAFT™ Suction Pumps if the packaging is broken.

**WARNING** DUTY CYCLE: The pump should be operated for 15 minutes ON with a 15-minute period OFF prior to re-operation.

**WARNING** Do not use in an area where flammable gases are present.

**WARNING** Regular periodic maintenance of the device is recommended annually.

**WARNING** Precautions for EMC safety should be observed. The device complies with EN60601-1-2:2015 for use in a Professional healthcare facility, however:

**WARNING** Electronic equipment in the vicinity of the device may affect its operation and potentially cause unpredictable operation of the device.

**WARNING** Wherever possible the device should be distanced from surrounding electromagnetic equipment and cables to this equipment in order to reduce possible electromagnetic interference.

**WARNING** The CRAFT™ Pump device power cable should only be connected to a correctly wired receptacle in order to avoid the risk of electrical shock and ONLY use the cable supplied by Rocket Medical.

**WARNING** When positioning the unit, ensure that access is available to the AC power cable located on the rear of the unit.

**WARNING** Fluids should not be allowed to enter the device as this may result in damage to the system.

**WARNING** User should be aware of the status of unit at all times during the procedure.

**WARNING** Only qualified personnel should service the device and the device must not be opened except by these personnel due to the risk of hazardous electrical shock and premature damage to the device. All service requirements should be referred to a Rocket Medical authorised representative.

*In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at [pncf@rocketmedical.com](mailto:pncf@rocketmedical.com) as well as to the competent health authority in the country that the user/patient resides.*

## Suction Pumps

## INSTRUCTIONS FOR USE

- WARNING** All equipment should be thoroughly cleaned after each use (Refer to Section 7).
- WARNING** Do not modify this equipment without written authorisation of the manufacturer.
- WARNING** The use of high suction levels may lead excessive fluid *flow rates* which may result in damage to the oocyte and reduced fertilisation rates. Damage to oocytes in harvesting systems is principally caused by *turbulent flow* which can lead to physical shearing stresses on the cumulus sufficient to denude or damage the fragile zona (Reeves et al 1989). Flow rate is a function needle set configuration and the suction applied.
- WARNING** Therefore, the -440mmHg (**Max**) suction setting must *only* be used to clear blockages from a needle set and must **NOT** be used in contact with the patient. Refer to Section 5 for recommended suction setting for given needle set configurations.

**2.6 BENEFIT RISK:** Rocket Medical plc has taken all necessary steps to ensure that residual risks associated with the use of CRAFT™ and CRAFT™ DUO-VAC™ Suction Pumps are reduced as far as possible through application of existing state of the art techniques in the design and manufacture of these medical devices to ensure safe usage. Rocket Medical plc concludes that the overall medical benefits of the CRAFT™ and CRAFT™ DUO-VAC™ Suction Pumps outweigh the possible risks when used according to the intended use.

**2.7 PERFORMANCE CLAIMS AND CLINICAL BENEFIT:** Aspiration of follicular fluid, oocytes and ovarian fluid and viable oocyte retrieval

**2.8 UNDESIRABLE SIDE EFFECTS:** The user should convey to the patient the potential risks of bleeding, infection and pain

**2.9 DESCRIPTION OF THE DEVICE:** The CRAFT™ and CRAFT™ DUO-VAC™ Suction Pumps have been developed to provide smooth, low volume vacuum at a pre-determined negative pressure. Vacuum is activated by a foot-operated toggle air switch controlled by the healthcare professional performing the gynaecological procedure. The range of vacuum is infinitely variable from 0 to -200 mmHg in medium vacuum mode and at a pre-set -440 mmHg in high vacuum mode.

**2.10 USER:** The intended user is a healthcare professional proficient in the undertaking of oocyte retrieval, working in accordance with local and national guidelines, or a trainee under the supervision of such a person.


**2.11 PATIENT POPULATION:** The intended patient demographics are adult patients undergoing oocyte aspiration or other gynaecological patients such as those undergoing oocyte retrieval for cryopreservation or those with ovarian hyperstimulation syndrome.

**2.12 ENVIRONMENT:** The CRAFT™ and CRAFT™ DUO-VAC™ Suction Pumps must be placed on a secure, level surface, away from sources of heat, water splashes, mists, or cooling vents. Do not expose to flammable gases.

	Temperature (°C)	Relative Humidity (%)
<b>Operating Limits</b>	+5 to +35	15 to 93
<b>Storage/Transport Limit</b>	-10 to +50	20 to 95

<b>WARNING</b> 	The device must be transported in temperatures +5°C to +35°C.
<b>WARNING</b> 	The device must be operated at relative humidity levels between 15% to 93%.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at [pncf@rocketmedical.com](mailto:pncf@rocketmedical.com) as well as to the competent health authority in the country that the user/patient resides.

<p>WARNING</p> 	<p>The device must be operated at ambient pressure levels between 70kPa to 106kPa.</p>
<p>WARNING</p> 	<p>Protect the device from ingress of liquid. Should any liquid enter the device, discontinue use immediately and refer to an authorised service agent.</p>
<p><b>Altitude</b></p>	<p>This device is intended for use below 2000 meters.</p>
<p>WARNING</p> 	<p>Device can cause explosion in the presence of flammable gases</p>
<p>WARNING</p> 	<p>The device is FRAGILE and must be transported in its original packaging to ensure protection. If the original packaging is not available please contact your local Customer Services Agent who will provide replacement packaging.</p>

**For Multiple Patient Multiple Use:** Intended to be used multiple times by multiple patients. Please follow cleaning instructions between use (See Section XX). At the end of each clinical session, turn off the device at the front panel and disconnect the device from mains power supply.

**Non-Sterile:** The devices is supplied non-sterile. It is not intended to be sterilised prior to use. Do not sterilise as doing so may compromise the structural integrity of the device, leading to device failure.

**MRI Safety:** MRI compatibility has not been verified.

## 2. GENERAL INFORMATION

### 3.1 SUPPLY VOLTAGE SELECTION

The device operates at a voltage 220-240 VAC @ 50Hz. 40VA or 110VAC @ 60Hz as appropriate to the model.

Ensure that the correct power cord is connected.

<p>WARNING</p> 	<p>Electric Shock hazard: The equipment is to be used only with electrical systems complying with all IEC, CEC and NEC requirements</p> <p>Do not immerse the device</p>
--	--

### 3.2 ELECTROMAGNETIC COMPATIBILITY

CRAFT™ Oocyte Aspiration Pumps comply with the electromagnetic compatibility (EMC) limits for medical devices as

*In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at [pncf@rocketmedical.com](mailto:pncf@rocketmedical.com) as well as to the competent health authority in the country that the user/patient resides.*

specified by BS EN60601-1-2:2015.

The device must be operated according to the instructions contained in this manual to ensure continued electromagnetic compatibility.

### 3.4 ESSENTIAL CONSUMABLES

R57685 Water Trap Set for CRAFT™ Suction Pumps

Glass water trap bottle, 2.0m patient connection tube with male luer connector. Supplied sterile, for single use in cartons of 10 units.

<p>WARNING</p> 	<p>This device must only be operated with approved consumable which are specifically designed to provide the correct flow rates and vacuum characteristics</p>
--	--

### 3.5 PACKAGING

The packaging has been carefully designed to allow secure transportation of the pump and its accessories. If the packaging is damaged, unintentionally opened before use or exposed to adverse environmental conditions check the device thoroughly before use.

### 3.6 SUCTION RANGE

Medium vacuum: 0-200mmHg<sup>-1</sup> in 20mmHg<sup>-1</sup> increments

High vacuum: 440mmHg<sup>-1</sup>

Accurate to ±5% full scale



### 3.7 CALIBRATION

CRAFT™ Pumps typically require little routine maintenance; however, they must be serviced and calibrated annually at a **Rocket Medical plc** approved service facility (see Section 8.2 Service Agents)


## 4. SYMBOLS AND SAFETY SIGNS USED ON THE DEVICE, LABELLING AND ACCESSORIES

Symbol or Safety Sign	Meaning	Symbol or Safety Sign	Meaning
	Manufacturer		Keep away from sunlight
	Authorized representative in the European Community/European Union		Keep dry
	Caution: Federal law restricts this device to sale by or on the order of a physician.		Temperature limit

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at [pncf@rocketmedical.com](mailto:pncf@rocketmedical.com) as well as to the competent health authority in the country that the user/patient resides.

	Date of manufacture		Consult the instructions for use
	Use-by date		Caution / Warning
<b>LOT</b>	Batch code	<b>SN</b>	Serial number
<b>REF</b>	Catalogue number		Does not contain natural rubber latex
<b>MD</b>	Medical device		Non-pyrogenic
	Do not use if package is damaged and consult instructions for use		Non-sterile
	Distributor		Importer
	Mains Power ON/OFF		This device is Type B
	Footswitch connection port		Dispose of this device in accordance with WEEE directive (2002/96/EC)
	Earthed		Read the manual
	Biological Risk		Usage Duty cycle: 15 minutes ON followed by 15 minutes OFF
	Humidity		Pressure
	Explosive Symbol		Electric Shock
<b>STERILE EO</b>	Sterilized using ethylene oxide	<b>PHT</b> <b>DEHP</b>	Contains or presence of phthalate

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at [pnrf@rocketmedical.com](mailto:pnrf@rocketmedical.com) as well as to the competent health authority in the country that the user/patient resides.

	Do not re-use		Single sterile barrier system
	Single sterile barrier system		Single sterile barrier system with protective packaging outside

## 5. OPERATING INSTRUCTIONS

CHECK: Do not use if the device is damaged.

CHECK: Ensure you have R57685 Water Trap Set for CRAFT™ Suction Pumps

CHECK: Ensure you have all consumables required for Oocyte Retrieval


PREPARE: Ensure the pump has been cleaned between patient use.

Ensure the Water Trap Set has not been used previously

### PROCEDURE:

#### CRAFT™ Suction Pumps (Ref: R29654 / R29655)

- I. (a) Unwrap the footswitch and attach the air tubing to the ports on the front of the pump casing. The footswitch normally remains attached to the unit in daily use unless removal for storage is required.

#### CRAFT™ DUO-VAC™ Suction Pumps (Ref: R29660 / R29661)

- (b) Unwrap the footswitch and attach the air tubing to the male/female ports on the front of the pump casing. The footswitch normally remains attached to the unit in daily use unless removal for storage is required.

*In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at [pncf@rocketmedical.com](mailto:pncf@rocketmedical.com) as well as to the competent health authority in the country that the user/patient resides.*

<b>WARNING</b> 	<b>Make sure the hose of the footswitch is not kinked</b>
--------------------	---

2. Connect the correct mains lead to an electrical supply 110-240VAC 50-60Hz. Model dependent.
3. Turn the power switch 0-I (front panel) to on. The green light illuminate.
4. Using aseptic technique, remove the R57685 Water Trap Set for CRAFT™ Suction Pumps from its packaging and attach the short tubing length to the water trap connection point of the CRAFT™ Suction Pump.

<b>WARNING</b> 	<b>This device must only be operated with approved consumable</b>
<b>WARNING</b> 	<b>Protect the device from ingress of liquid. Should any liquid enter the device, discontinue use immediately and refer to an authorised service agent</b>
<b>WARNING</b> 	<b>The R57685 Water Trap Set for CRAFT™ Suction Pumps is for single-use and must be replaced for each patient.</b>

5. Pass the longer patient connection tube into the operative field.
6. Refer to Table I for recommended vacuum values for various needle and tube set combinations.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at [pncf@rocketmedical.com](mailto:pncf@rocketmedical.com) as well as to the competent health authority in the country that the user/patient resides.

	Single Lumen		Double Lumen	
	16G	17G	16G	17G
Tube Set Length	Recommended Vacuum (-mmHg)			
55cm	80	110	130	150
70cm	90	130	150	170
90cm	100	150	170	190

Table I Recommended vacuum settings

*These values are provided for the clinician's consideration; for guidance purposes only.*

<b>WARNING</b> 	<b>Ensure the set vacuum level is appropriate to the patient's needs taking into consideration the needle and tube set configuration.</b>
---	---

7. Occlude the patient tube set, distal to the water trap.
8. Activate the footswitch.
9. Rotate the control knob clockwise to increase the vacuum level to the desired value on the gauge.
10. To decrease the set value, rotate anticlockwise.
11. When the desired value is reached, release the footswitch and tubing set occlusion.
12. Connect the patient filter set to a suitable luer fitting (indicated) such as the Oxford Tube Set found on Rocket Medical Embryon® Single or Double Lumen Oocyte Aspiration Needles (not supplied). Attach a collection tube to the bung (not supplied)
13. Proceed to follow the directions provided in the manufacturer's IFU of the aspiration needle for the oocyte retrieval procedure.

## 6. TROUBLESHOOTING

### 6.1 BLOCKED VACUUM LINES

If vacuum lines become blocked due to debris or viscous fluids, a temporary high vacuum (-440mmHg) can be obtained to unblock the needle and tubing set by following the procedures below.

<b>WARNING</b> 	<b>HIGH VACUUM MODE</b> must NOT be used to aspirate oocytes as this may result in damage to the oocyte and lead to reduced fertilisation rates.
---	---

### CRAFT™ DUO-VAC™

To activate the HIGH vacuum, depress the WHITE HIGH VACUUM Pedal – the pump will immediately deliver a vacuum of -440mmHg.

*In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at [pncf@rocketmedical.com](mailto:pncf@rocketmedical.com) as well as to the competent health authority in the country that the user/patient resides.*

### CRAFT™

To activate the HIGH vacuum, occlude the RED button and depress the footswitch – the pump will immediately deliver a vacuum of -440mmHg. Removing the occlusion will cause the vacuum to revert to its previously set level.

### 6.2 WATER TRAP SET

If water, media or any other material is allowed to enter the water trap chamber, the pump must be stopped **immediately**. Replace with a R57685 Water Trap Set for CRAFT™ Suction Pumps. The pump must NOT be operated if fluid is present in the water trap. Should the water trap become full and there is a risk that fluid has entered the pump, it must be withdrawn from service and returned to an authorised service agent for inspection and repair

<p><b>WARNING</b></p> 	<p><b>This device must only be operated with approved consumable which are specifically designed to provide the correct flow rates and vacuum characteristics</b></p> <p><b>Should any liquid enter the device, discontinue use immediately and refer to an authorised service agent</b></p>
---	--

## 7. DEVICE MAINTENANCE

At the end of each clinical session, turn off the device at the front panel and disconnect the device from mains power supply. Using an aqueous solution of 70% alcohol (IMS or isopropyl BP), moisten a cloth and wipe all external surfaces of the device. If the surface has become contaminated with proteinaceous material, remove with a light detergent solution before surfacing cleaning with alcohol.

*In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at [pncf@rocketmedical.com](mailto:pncf@rocketmedical.com) as well as to the competent health authority in the country that the user/patient resides.*

Do NOT use a 100% alcohol or any other solvent to clean the device as this may cause damage to the casing surface and display.

Prevent any fluid from entering the device.

<p><b>WARNING</b></p> 	<p><b>Electric Shock hazard: Do not immerse the device.</b></p>
<p><b>WARNING</b></p> 	<p><b>Do NOT attempt to sterilise the device</b></p>
<p><b>WARNING</b></p> 	<p><b>Protect the device from ingress of liquid. Should any liquid enter the device, discontinue use immediately and refer to an authorised service agent</b></p>

## 8. SERVICE INFORMATION

### 8.1 RETURNING THE PUMP

**All devices to be returned must be prepared as described below for the protection of the servicing team and for safety during transport.**

Surface clean the pump as described in the Section above.

Seal in a plastic bag and seal within a second plastic bag.

Place in the original packaging.

Enclose the following information:

Contact name

Centre address

Decontamination Certificate

Description of the fault or service required

Accompanying Order to authorise servicing – contact your local Customer Services Team for details.

CRAFT™ Pumps typically require little routine maintenance; however, they must be serviced and calibrated annually at a

**Rocket Medical plc** approved service facility.

Failure to service the pump at the indicated intervals may invalidate the Warranty.

UK & European Service Agents:

**\*IMPORTANT:**

<p><b>WARNING</b></p> 	<p><b>No user serviceable parts inside.</b></p>
<p><b>WARNING</b></p> 	<p><b>Do NOT include used consumables as these pose a significant contamination risk</b></p>
<p><b>WARNING</b></p> 	<p><b>A decontamination certificated MUST be included with every returned pump. Repair or servicing cannot be commenced unless the service agent is in possession of this certificate.</b></p>

### 8.2 SERVICE AGENTS

*In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at [pncf@rocketmedical.com](mailto:pncf@rocketmedical.com) as well as to the competent health authority in the country that the user/patient resides.*

## Suction Pumps

## INSTRUCTIONS FOR USE

The CRAFT™ and CRAFT™ DUO-VAC™ Suction Pumps typically require little routine maintenance. However, they must be serviced and calibrated **annually** at a service facility approved by Rocket Medical plc. Failure to service the pump at the indicated intervals may invalidate the Warranty.

**UK & Ireland:****IVF Synergy Ltd.**

1st Floor Suite  
Maritime House  
Discovery Quay  
Falmouth  
TR11 3XA  
UNITED KINGDOM  
Tel: +44 (0) 1872 487224  
Email: [service@ivfsynergy.com](mailto:service@ivfsynergy.com)

**Australian Service Agent:****Rocket Medical Pty Ltd**

Suite 6,  
157 Gordon Street,  
Port Macquarie NSW 2444  
AUSTRALIA  
Email:  
[aus@rocketmedical.com](mailto:aus@rocketmedical.com)

**Rest of World:****Charles Austen Pumps**

Charles Austen Pumps Ltd,  
Royston Road,  
Byfleet, KT14 7NY.  
UNITED KINGDOM  
Tel: +44 (0) 1932 355277  
Email:  
[Sales@charlesausten.com](mailto:Sales@charlesausten.com)

**UK Customer Service:**

Rocket Medical plc.  
Sedling Road  
Washington  
NE38 9BZ  
UNITED KINGDOM  
Tel: +44 (0) 191 419 6988  
Fax: +44 (0) 191 419 6989  
Email:  
[customerservices@rocketmedical.com](mailto:customerservices@rocketmedical.com)

**US Office:****Rocket Medical**

50 Corporate Park Drive. Suite  
890.  
Pembroke. MA. 02359. USA  
Tel: +1 781 749 6223  
Email:  
[usa@rocketmedical.com](mailto:usa@rocketmedical.com)

## 9. WARRANTY

CRAFT™ Oocyte Aspiration Pumps are sold by Rocket Medical plc under the warranties set forth in the following paragraphs. Such warranties are extended only with respect to the purchase of the Products directly from Rocket Medical plc as new merchandise and are extended to the first Buyer thereof, other than for resale.

For a period of TWENTY-FOUR (24) months from the date of shipment the Products are warranted to be free from functional defects in materials and workmanship and to conform to the description of the Products contained in the operating manual and accompanying labels, provided the same is properly operated under conditions of normal use, that annual maintenance and service is performed at an authorised Rocket Medical plc service facility.

Removal of any QC seal voids the warranty.

The foregoing warranties shall not apply if the Products have been repaired other than by Rocket Medical plc or other than in accordance with written instructions provided by Rocket Medical plc, or altered by anyone other than Rocket Medical plc, or if the Products have been subject to misuse, negligence, or accident.

Rocket Medical plc's sole and exclusive obligation and Buyer's sole and exclusive remedy under the above warranties is limited to repairing or replacing, free of charge, at Rocket Medical plc's option, Products, which are reported to Rocket Medical plc by mail, telephone or email and which, if so advised by Rocket Medical plc, is thereafter returned with a statement of the observed deficiency, not later than seven (7) days after the expiration date of the warranty, to Rocket Medical plc during normal business address, transport charges prepaid and which, upon Rocket Medical plc's examination, is not found to conform with the above warranties.

Rocket Medical plc shall not be otherwise liable for any damages including but not limited to incidental damages, consequential damages or special damages.

There are no express or implied warranties which extend beyond the warranties herein above set forth. Rocket Medical plc makes no warranty of merchantability or fitness for a particular purpose with respect to the Products or parts thereof.

*In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at [pncf@rocketmedical.com](mailto:pncf@rocketmedical.com) as well as to the competent health authority in the country that the user/patient resides.*

## 10. TECHNICAL SPECIFICATIONS

### 10.1 CLASSIFICATION

IEC 60601-1

Type of protection against electric shock: Class I. Degree of protection against electric shock: Type B. Vacuum type: high vacuum/low volume.

Suitable for continuous operation.

Not suitable for use in the presence of flammable gases.

Not suitable for use in conditions which expose the device to the ingress of water. Not suitable for sterilisation.


### 10.2 SPECIFICATIONS

Power Input to Pump: 220-240 VAC @ 50Hz. 40VA Maximum current: 2.5A

Consumption: 0.6KW/h Environmental conditions:

- Temperature +5°C to +35°C
- Atmospheric Pressure Range: 700-1060hPA Dimensions:

:

	<p><b>Dimensions:</b></p> <ul style="list-style-type: none"> <li>• W - 264mm</li> <li>• H - 124mm</li> <li>• D - 164mm</li> </ul>	<p><b>Weight:</b></p> <ul style="list-style-type: none"> <li>• Unit – 1.9Kg</li> <li>• Foot Switch: <ul style="list-style-type: none"> <li>○ 0.57Kg – DUO-VAC™ (plastic)</li> <li>○ 0.34kg – Craft Suction (plastic)</li> <li>○ 0.51Kg – DUO-VAC™ (metal)</li> <li>○ 0.38kg – Craft Suction (metal)</li> </ul> </li> </ul>
---	---	--

- Unit – 1.9Kg
  - Medium vacuum: 0 to -200mmHg in -20mmHg increments
  - High vacuum: -440mmHg
- Accurate to ±5% full scale

### 10.3 EMC TABLES:

**Table 1: Guidance and Manufacturer's Declaration – Electromagnetic Emissions**

Emissions Test	Compliance	Electromagnetic environment - guidance
RF Emissions CISPR 11	Group 1	The CRAFT™ Aspiration Pumps generate RF signals for its internal function. Therefore, its RF emissions are very low and not likely to cause any interference in nearby electronic equipment.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at [pncf@rocketmedical.com](mailto:pncf@rocketmedical.com) as well as to the competent health authority in the country that the user/patient resides.

RF CISPR 11	Class A	CRAFT™ Aspiration Pumps are suitable for use in all establishments, other than domestic and those connected directly to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

**Table 2: Guidance and Manufacturer's Declaration - Electromagnetic Immunity**

CRAFT™ Aspiration Pumps (the device) is intended for use within the electromagnetic environment specified below. The operator of a device should ensure that it is used within such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic environment Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+/- 6kV contact +/- 8kV air	+/- 6kV contact +/- 8kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient burst IEC 61000-4-4	+/- 2kV for power supply lines	+/- 2kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	+/- 1kV differential mode +/- 2kV common mode	+/- 1kV differential mode +/- 2kV common mode	Mains power quality should be that of a typical commercial or hospital environment
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at [pncf@rocketmedical.com](mailto:pncf@rocketmedical.com) as well as to the competent health authority in the country that the user/patient resides.

Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-11	<5% Ut (>95% dip in Ut) for 0.5 cycle 40% Ut (60% dip in Ut) for 5 cycles 70% Ut (30% dip in Ut) for 25 cycles <5% UT (>95% dip in Ut) for 5 seconds	<5% Ut (>95% dip in Ut) for 0.5 cycles 40% Ut (60% dip in Ut) for 5 cycles 70% Ut (30% dip in Ut) for 25 cycles <5% UT (>95% dip in Ut) for 5 seconds	Mains power quality should be that a typical commercial or hospital environment. If the operator of the device requires continued operation during mains interruptions, it is recommended that the device be powered from an uninterruptable power
Note: Ut is the a.c. mains voltage prior to application of the test level.			

**Table 3: Guidance and Manufacturer’s Declaration - Electromagnetic Immunity**

The device is intended for use within the electromagnetic environment specified below. The customer or the operator of a should ensure that it is used within such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic environment Guidance
			Portable and mobile RF communications equipment should be no closer to any part of the PSU device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter
			Recommended Separation Distance
Conducted Rf	3 Vrms	3Vrms	$D=1.2\sqrt{P}$
IEC 61000-4-6	150 kHz to 80MHz		
Radiated RF	3V/m	3V/m	$D=1.2\sqrt{P}$ 80 MHz to 800 MHz
IEC 61000-4-3	80 MHz to 2.5 GHz		$D=1.2\sqrt{P}$ 800 MHz to 300 GHz
			Where P is the maximum power output [power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m) Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey a, should be less than the compliance level in each frequency in each frequency range.
			Interference may occur in the vicinity of equipment marked with the following symbol

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at [pnrf@rocketmedical.com](mailto:pnrf@rocketmedical.com) as well as to the competent health authority in the country that the user/patient resides.

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and /or people.

a. Field strengths from fixed transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.

b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

**Table 4: Recommended Separation Distances Between Portable and Mobile RF Communications Equipment and the Oocyte Aspiration Pump.**

The device is intended for use in the electromagnetic environment specified below.

The operator of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Read Maximum Output Power of Transmitter	Separation Distance According to Frequency of Transmitter m		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
W	$d = 1.2 \sqrt{P}$	$d = 1.2 \sqrt{P}$	$d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	3.8	3.8	7.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distances d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and/or people.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at [pncf@rocketmedical.com](mailto:pncf@rocketmedical.com) as well as to the competent health authority in the country that the user/patient resides.

## II. DISPOSAL

This device, its accessories and the consumables used with it, should be handled and disposed of in accordance with policy of the healthcare setting and with regard to all applicable regulations, including but without limitation to, those pertaining to human health & safety and care of the environment. Particular care to be taken with sharp objects such as needles due to the risk of infection, cuts and environmental hazards.

**Waste Electrical and Electronic Equipment:** In Europe, disposal of this device must be undertaken with regard to the WEEE directive (2002/96/EC). Do not dispose of this device in municipal waste.

*In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at [pncf@rocketmedical.com](mailto:pncf@rocketmedical.com) as well as to the competent health authority in the country that the user/patient resides.*

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