
Product Name: OOCYTE ASPIRATION PUMP

ZDOCK193 Rev. 34

Section	Language	
1	en	Oocyte Aspiration Pump Instructions For Use
2	cn	卵母细胞抽吸泵 使用说明
3	de	Oozytenaspirationspumpe Gebrauchsanweisung
4	dk	Oocyte Aspirationspumpe Brugsanvisning
5	es	Bomba De Aspiración De Ovocitos Instrucciones De Uso
6	gr	ΑΝΤΛΙΑ ΑΝΑΡΡΟΦΗΣΗΣ ΩΟΚΥΤΤΑΡΩΝ ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
7	nl	OÖCYT ASPIRATIEPOMP GEBRUIKSAANWIJZING
8	pl	Pompa Aspiracyjna Do Oocytów Instrukcja Obsługi
9	se	Pump För Aspiration Av Oocyter Bruksanvisning
10	vi	Hướng dẫn sử dụng bơm hút noãn bào
11	it	Istruzioni per l'uso Pompa di aspirazione degli ovociti
12	cz	Aspirační pumpa oocytů

ROCKET MEDICAL PLC
Sedling Road, Washington,
NE38 9BZ, UK

ROCKET MEDICAL GmbH
Am Rosengarten 48,
15566 Schöneiche.
Germany

ROCKET MEDICAL GmbH
Am Rosengarten 48,
15566 Schöneiche.
Germany

ROCKET MEDICAL LLC
50 Corporate Park Drive #890.
East Pembroke. MA. 02359.
USA

Australian Sponsor:

Rocket Medical Pty Ltd
Suite 6, 157 Gordon Street,
Port Macquarie.
NSW 2444. Australia

COPYRIGHT

This manual contains information that is subject to copyright. All rights reserved.

This manual should not be photocopied, duplicated or distributed completely or in part without the approval of Rocket Medical plc

Oocyte Aspiration Pump

CONTENTS

1. SAFETY INSTRUCTIONS	2
2. ABOUT OOCYTE ASPIRATION PUMP	3
2.1 INTENDED USE	4
2.2 INDICATIONS	4
2.3 CONTRAINDICATIONS	4
2.4 WARNINGS:	4
2.5 BENEFIT RISK	5
2.6 PERFORMANCE CLAIMS AND CLINICAL BENEFIT	5
2.7 UNDESIRABLE SIDE EFFECTS	5
2.8 DESCRIPTION OF THE DEVICE:	5
2.9 USER	5
2.10 PATIENT POPULATION	5
2.11 ENVIRONMENT	5
3. GENERAL INFORMATION	6
3.1 TERMINOLOGY	6
3.2 POWER SUPPLY	6
3.3 SUPPLY VOLTAGE SELECTION	6
3.4 ELECTROMAGNETIC COMPATIBILITY	7
3.5 ESSENTIAL CONSUMABLES	7
3.6 PACKAGING	7
3.7 SUCTION RANGE	7
4. SYMBOLS AND SAFETY SIGNS USED ON THE DEVICE, ACCESSORIES, PSU AND LABELLING	8
5. OPERATING INSTRUCTIONS	9
6. PUMP SET-UP PROCEDURE SCHEMATIC	12
7. DEVICE MAINTENANCE	12
8. SERVICE INFORMATION	13
8.1 SERVICE INDICATOR	13
8.2 RETURNING THE PUMP FOR SERVICE	13
8.3 SERVICE AGENTS	13
9. WARRANTY	14
10. TECHNICAL SPECIFICATIONS	14
10.1 CLASSIFICATION	14
10.2 SPECIFICATIONS	15
10.3 EMC TABLES:	15
11. DISPOSAL	18

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

I. SAFETY INSTRUCTIONS

This manual describes the operation and intended use of the device and its associated consumables. It is essential that this document is used to familiarise yourself with the correct function and operation of the device before use.

Failure to follow these instructions may result in serious injury to the patient or operator and can lead to damage or breakdown of the device. In case the device fails during an operation, a replacement device and replacement disposables should be kept within reach so that the operation may be completed.

This manual does not provide a detailed description of the oocyte harvesting procedure and is not intended as a training guide for users inexperienced in the technique.

The device must be used with R57686 Filter Set for Oocyte Aspiration Pump. The usage of non-approved tubing or filter sets may impair pump performance, lead to increased risk to patients and harvested oocytes, and will invalidate the warranty.

The R57686 Filter Set for Oocyte Aspiration Pump is designed to prevent fluid contamination of the vacuum pump.

If the suction pump has been used with a non-approved filter set or there is any evidence or suspicion that the pump may have been contaminated with fluid during use, it must be removed from service and returned for examination immediately. Please contact your nearest Service Centre for advice.

Refer all servicing to the manufacturer's authorised service agent.

WARNING 	READ THIS MANUAL CAREFULLY: Please familiarise yourself with the contents of this manual before attempting to use the device. Failure to observe these instructions may result in damage to the pump or cause injury to the patient or user. This device should only be used by suitably qualified personnel.
WARNING 	Electric Shock hazard: The equipment is to be used only with electrical systems complying with all IEC, CEC and NEC requirements
WARNING 	Any adjustment, modification or repairs to the equipment should be carried out by authorised service agents. This device should only be operated by appropriately qualified personnel Ensure the set vacuum level is appropriate. No user serviceable parts inside
DISPOSAL 	Disposal of this device must be undertaken with regard to the WEEE directive (2002/96/EC).
WARNING 	This device must only be operated with approved consumables
WARNING 	Protect the device from ingress of liquid. Should any liquid enter the device, discontinue use immediately and refer to an authorised service agent
WARNING 	Device can cause explosion in the presence of flammable gases

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

2. ABOUT OOCYTE ASPIRATION PUMP

Product Names: Oocyte Aspiration Pump

Product Codes: R29700

Device Image:

Figure 1 Oocyte Aspiration Pump

1. O/I 12V Power On/Off switch (not visible)
2. Suction Display in -mmHg
3. Patient connection port – only for use with R57686 Filter Set for Oocyte Aspiration Pump
4. Footswitch connection ports (not visible)
5. LEDs (Top to Bottom)
 - Pre-set (Max) Suction Indicator LED (Orange, -500mmHg)
 - User Set Suction Indicator LED (Blue, -50 to -300mmHg)
 - Power-On Indicator LED (Green, 12VDC)
6. Touch sensitive Suction Control Dial – clockwise to increase, anticlockwise to decrease the set value
7. Service Indicator LED (Yellow)
8. Dual footswitch – air controlle

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

2.1 INTENDED USE: Intended for the generation of low volume/high (vacuum) suction of between -50 to -300 mmHg to permit the aspiration of follicular fluid, oocytes and ovarian fluid.

2.2 INDICATIONS: Indicated for use as part of the treatment for infertility relating to IVF and other related procedures. This device may only be used by, or under the supervision of, appropriately trained personnel in conjunction with clinical practice guidelines.

2.3 CONTRAINDICATIONS:

The patient should be made aware of:

- Not intended for use where ovarian aspiration or the aspiration of ovarian fluid is contraindicated:
- Patients with vaginal infections and or sexually transmitted diseases as this increases the risk of transmission of infection from the pump to the next patient
- Suspected pregnancy as this could result in miscarriage or fetal injury
- Where the pelvic anatomy cannot be accurately assessed by ultrasonography as this could result in damage to other organs and vessels
- Not intended for surgical suction or other applications as the diameter of the tubing is too narrow for adequate drainage.

2.4 WARNINGS:

WARNING	Users should be familiar with and adhere to all warnings, cautions and instructions for use that are labelled on the device and included in the User Manual.
WARNING	This device should only be used by, or under the supervision of, appropriately trained personnel and in conjunction with current local clinical practice guidelines.
WARNING	The device is not approved for connection to any other device except the R57686 Filter Set for Oocyte Aspiration Pump
WARNING	The R57686 Filter set for Oocyte Aspiration Pump is single use and its instructions for use must be followed at all times.
WARNING	Do NOT use the R57686 Filter set for Oocyte Aspiration Pump if the packaging is broken.
WARNING	Do not use in an area where flammable gases are present.
WARNING	Regular periodic maintenance of the device is recommended as indicated by the service indicator.
WARNING	Precautions for EMC safety should be observed. The device complies with EN60601-1-2:2015 for use in a Professional healthcare facility however:
WARNING	Electronic equipment in the vicinity of the device may affect its operation and potentially cause unpredictable operation of the device.
WARNING	Wherever possible the device should be distanced from surrounding electromagnetic equipment and cables to this equipment in order to reduce possible electromagnetic interference.
WARNING	The Oocyte device power cable should only be connected to a correctly wired receptacle in order to avoid the risk of electrical shock and ONLY use the cable supplied by Rocket Medical.
WARNING	When positioning the unit, ensure that access is available to the DC power cable located on the rear of the unit.
WARNING	To avoid damage to the unit only use the approved mains power adaptor supplied with the unit, Model Number: MPU30.
WARNING	Fluids should not be allowed to enter the device as this may result in damage to the system.
WARNING	User should be aware of the status of unit at all times during the procedure.
WARNING	Only qualified personnel should service the device and the device must not be opened except by these personnel due to the risk of hazardous electrical shock and premature damage to the device. All service requirements should be referred to a Rocket Medical authorised representative.
WARNING	All equipment should be thoroughly cleaned after each use (refer to section 6 Cleaning Guidelines).
WARNING	Do not modify this equipment without written authorisation of the manufacturer.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

WARNING	The use of high suction levels may lead excessive fluid <i>flow rates</i> which may result in damage to the oocyte and reduced fertilisation rates. Damage to oocytes in harvesting systems is principally caused by <i>turbulent flow</i> which can lead to physical shearing stresses on the cumulus sufficient to denude or damage the fragile zona (Reeves et al 1989). Flow rate is a function needle set configuration and the suction applied. Therefore, the -500mmHg (Max) suction setting must only be used to clear blockages from a needle set and must NOT be used in contact with the patient. Refer to Section 6 for recommended suction setting for given needle set configurations.
---------	---

2.5 BENEFIT RISK: Rocket Medical plc has taken all necessary steps to ensure that residual risks associated with the use of Oocyte Aspiration Pump are reduced as far as possible through application of existing state of the art techniques in the design and manufacture of these medical devices to ensure safe usage. Rocket Medical plc concludes that the overall medical benefits of the Oocyte Aspiration Pump outweigh the possible risks when used according to the intended use.

2.6 PERFORMANCE CLAIMS AND CLINICAL BENEFIT: Aspiration of follicular fluid, oocytes and ovarian fluid and viable oocyte retrieval

2.7 UNDESIRABLE SIDE EFFECTS: The user should convey to the patient the potential risks of bleeding, infection and pain

2.8 DESCRIPTION OF THE DEVICE: The Oocyte Aspiration Pump has been developed to provide smooth, low volume/high suction (vacuum) at a pre-determined negative pressure. Suction is activated by a foot operated toggle air switch controlled by the surgeon performing the oocyte collection. The range of suction is variable from – 50 to - 300mmHg and at a pre-set -500mmHg in 'Max' suction mode. The Oocyte Aspiration Pump requires a disposable filter set.

2.9 USER: The intended user is a healthcare professional proficient in the undertaking of oocyte retrieval, working in accordance with local and national guidelines, or a trainee under the supervision of such a person.

2.10 PATIENT POPULATION: The intended patient demographics are adult patients undergoing oocyte aspiration or other gynaecological patients such as those undergoing oocyte retrieval for cryopreservation or those with ovarian hyperstimulation syndrome.

2.11 ENVIRONMENT: The Oocyte Aspiration Pump must be placed on a secure, level surface, away from sources of heat, water splashes, mists, or cooling vents. Do not expose to flammable gases.

	Temperature (°C)	Relative Humidity (%)
Operating Limits	+5 to +40	15 to 93
Storage/Transport Limit	-10 to +50	20 to 95

WARNING 	The device must be transported in temperatures +5°C to +40°C.
WARNING 	The device must be operated at relative humidity levels between 15% to 93%.
WARNING 	The device must be operated at ambient pressure levels between 70kPa to 106kPa.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

WARNING 	Protect the device from ingress of liquid. Should any liquid enter the device, discontinue use immediately and refer to an authorised service agent.
Altitude	This device is intended for use below 2000 meters.
WARNING 	Device can cause explosion in the presence of flammable gases
WARNING 	The device is FRAGILE and must be transported in its original packaging to ensure protection. If the original packaging is not available please contact your local Customer Services Agent who will provide replacement packaging.

For Multiple Patient Multiple Use: Intended to be used multiple times by multiple patients. Please follow cleaning instructions between use (See Section XX). At the end of each clinical session, turn off the device at the rear panel and disconnect the device from mains power supply.

Non-Sterile: The device is supplied non-sterile. It is not intended to be sterilised prior to use. Do not sterilise as doing so may compromise the structural integrity of the device, leading to device failure.

MRI Safety: MRI compatibility has not been verified.

3. GENERAL INFORMATION

3.1 TERMINOLOGY Throughout this manual the term “suction” is used to denote aspiration of fluid with vacuum or negative pressure. The use of the term ‘suction’ is in observance of the requirements of the ISO6061 -I Standard as amended.

3.2 POWER SUPPLY: The device is only for use with power supply Model Number: MPU30. Attachment of any other power supply may severely damage the device.

CAUTION 	Disconnection from the mains supply can only be achieved with the removal of the mains power lead from the wall socket.
--	--

3.3 SUPPLY VOLTAGE SELECTION The device operates at a voltage of 12VDC supplied via a universal medical grade power supply unit (PSU). Input voltage: 12VDC Universal power supply range 100V to 240 VAC @ 50 / 60Hz Ensure that the correct power cord is connected.

WARNING 	Electric Shock hazard: The equipment is to be used only with electrical systems complying with all IEC, CEC and NEC requirements Do not immerse the device
--	--

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

3.4 ELECTROMAGNETIC COMPATIBILITY The Oocyte Aspiration Pump complies with the electromagnetic compatibility (EMC) limits for medical devices as specified by EN60601-1-2:2015. These limits are designed to provide a reasonable degree of protection against harmful interference found in typical medical installations. Medical electrical equipment requires special precautions regarding EMC and the device must be installed, positioned and operated according to the instructions contained in this manual to ensure continued electromagnetic compatibility. It is possible that high levels of radiated or conducted radio-frequency electromagnetic interference (EMI) from portable and mobile RF communications equipment or other strong or nearby radio-frequency sources could result in performance disruption of the suction pump.

Evidence of EMC interference may include erratic digital displays, an inability to correctly set the desired suction, the device failing to operate, or other incorrect functioning.

If this occurs, stop using the aspiration pump and contact Customer Services or Rocket® Medical authorised distributor.

3.5 ESSENTIAL CONSUMABLES R57686 Filter Set for Oocyte Aspiration Pump Set 5µm filter set, 2.5m patient connection tube with male luer connector. Supplied sterile, for single use in cartons 10 units.

WARNING 	This device must only be operated with approved consumable which are specifically designed to provide the correct flow rates and <i>protect the device from ingress of liquid.</i>
WARNING 	Do not use Aspiration set if the packaging is broken.
WARNING 	Should any liquid enter the device, discontinue use immediately and refer to an authorised service agent

3.6 PACKAGING The packaging has been carefully designed to allow secure transportation of the pump and its accessories. If the packaging is damaged, unintentionally opened before use or exposed to adverse environmental conditions check the device thoroughly before use.

CAUTION 	After unpacking, re-assemble and retain the packaging for transport for servicing when required.
---	--

3.7 SUCTION RANGE

Suction Ranges:

User Set: -50mmHg to -300mmHg in 5mmHg increments

Pre-set: -500mmHg - nominal


Accurate to ±5% of full scale

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

4. SYMBOLS AND SAFETY SIGNS USED ON THE DEVICE, ACCESSORIES, PSU AND LABELLING

Symbol or Safety Sign	Meaning	Symbol or Safety Sign	Meaning
	Manufacturer		Keep away from sunlight
	Authorized representative in the European Community/European Union		Keep dry
	Caution: Federal law restricts this device to sale by or on the order of a physician.		Temperature limit
	Date of manufacture		Medical device
	Use-by date		Consult the instructions for use
	Batch code		Caution
	Catalogue number		
	Serial number		Does not contain natural rubber latex
	Non-sterile		Contains or presence of phthalate
	Do not use if package is damaged and consult instructions for use		Non-pyrogenic
	Distributor		Importer
	Mains Power ON/OFF		This device is Type B
	Footswitch connection port		Dispose of this device in accordance with WEEE directive (2002/96/EC)
	Earthed		Read the manual before connection and use

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

	Service interval indicator		Filter set connection point
	Suction display touch control symbol		Indicates activation of the pre-set (Max) -500 mmHg mode
	Indicates activation of the user-set -50 to -300 mmHg mode		Dispose of this device in accordance with WEEE directive (2002/96/EC)
	Electric Shock Risk		Device is only for use indoors
	Humidity		Pressure
	Class II Power Supply Unit		12V DC connection polarity
	PSU conforms to EN60601-1 & IEC60601-1		CE Mark
	Single sterile barrier system		Single sterile barrier system with protective packaging outside
	Sterilized using ethylene oxide		Do not resterilise
	Do not reuse		Biological Risk

5. OPERATING INSTRUCTIONS


CHECK: Do not use if the device is damaged.

CHECK: Ensure you have R57686 Water Filter Set for Oocyte Aspiration Pump

CHECK: Ensure you have all consumables required for Oocyte Retrieval

PREPARE: Ensure the pump has been cleaned between patient use.

Ensure the Filter Set packaging is not damaged or broken and has not been used on a previous patient.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

1. Unwrap the footswitch and attach the air tubing to the male/female ports on the underside of the pump casing. The footswitch is normally left attached to unit in daily use unless removal for storage is required. (Figure 2)
2. Connect the correct mains lead to the 12V power supply module - Model number: MPU30 and connect to an electrical supply 100V to 240 VAC @ 50 / 60Hz. The power supply module is self-regulating to generate 12VDC.
3. Turn the power switch 0-I ☐ (rear panel) to on. The green LED (Figure 3) on the front panel will illuminate. (Figure 3)
4. The main display will display an initial value of -100 mmHg.
5. Using aseptic technique, unpack the R57686 Filter Set for Oocyte Aspiration Pump and attach the short length of silicone tubing to the spigot (Figure 4). The tubing should be a light push-on fit.*
6. Pass the longer patient connection tube into the operative field. The filter set is for single-use only and must be replaced for each patient.

Figure 2

Figure 3

Figure 4

WARNING 	If the filter material becomes contaminated with fluid it will occlude thereby protecting the pump mechanism. Immediately replace the filter set if it becomes contaminated and between every patient.
WARNING 	The pump must NOT be operated without the correct filter set attached. The hydrophobic filter set is specially designed to prevent ingress of fluid to the pump.
WARNING 	This device must only be operated with approved consumables Operation without the correct filter set will invalidate the Warranty.
WARNING 	ELECTRIC SHOCK HAZARD The equipment is to be used only with electrical systems complying with all IEC, CEC and NEC requirements.
WARNING 	Make sure the footswitch hose is not kinked

7. On Start-Up, the pump is pre-set to -100 mmHg suction. To increase this value, use a circular motion, lightly drawing a finger clockwise over the outer zone of the control dial in the direction of the arrows. (Figure 5)
8. Refer to Table 1 below for recommend values for various needle and tube set combinations
9. The active dial zone will sense the motion and increase the display in -5mmHg increments.
10. To decrease the set value, use an anticlockwise circular motion.
11. When the desired value is reached, hold the finger still for 2 seconds when an audible 'beep' will confirm the new setting.
12. If the new setting is not correctly saved, the value will return to its previous set value in 10 seconds.
13. If the new value has not been stored, repeat from Step 7.

Figure 5

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

	Single Lumen		Double Lumen	
	16G	17G	16G	17G
Tube Set Length	Recommended Suction -mmHg			
55cm	80	110	130	150
70cm	90	130	150	170
90cm	100	150	170	190

WARNING 	Ensure the set suction level is appropriate to the patient's needs and has taken into consideration the needle and tube set configuration
---	---

14. Connect the patient filter set to a suitable luer fitting (indicated) such as that found on the Oxford Tube Set found on Single and Double Lumen Oocyte Aspiration Needles and attach a collection tube to the bung (Figure 6)
15. The luer connection is also suitable for connection to other standard ISO 6% luer taper adaptors and fittings.
16. The footswitch control operates on a toggle switch. Depress a pedal once for ON and depress again for OFF. A tone will sound to indicate activation of the footswitch.
17. To activate the suction pump to deliver the set value – depress the WHITE pedal once.
18. Once activated, the pump will commence to deliver the suction level previously set. The display indicates the actual suction set by the user.
19. The blue LED will illuminate and flash rapidly to confirm suction (-50 to -300 mmHg) is active. (Figure 8)
20. To stop suction, depress the pedal once again. The pump will immediately stop and vent to air and the user set suction LED will go out.
21. To activate the pre-set Max suction, depress the BLACK Pedal once – the pump will immediately deliver a nominal suction of -500 mmHg. (Figure 9)
22. Both suction indicator LEDs will illuminate.
23. To stop Max suction, depress the BLACK pedal once or depress the WHITE pedal to continue to evacuate at the previously set suction level.
24. Follow the manufacturer's instructions provided with the aspiration needle for the oocyte retrieval procedure.
25. If the pump is powered off between uses, please note that it reverts back to the pre-set value of -100mmHg when it is powered back on.


Figure 6


Figure 7

Figure 8


Figure 9

WARNING 	High Suction Mode must NOT be used to aspirate oocytes as this may result in damage to the oocyte and lead to reduced fertilisation rates. High Suction settings may only be used to clear blockages from a needle set and must NOT be used in contact with the patient.
---	--

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

6. PUMP SET-UP PROCEDURE SCHEMATIC

7. DEVICE MAINTENANCE

At the end of each clinical session, turn off the device at the rear panel and disconnect the device from mains power supply.

Using an aqueous solution of 70% alcohol (IMS or isopropyl BP), moisten a cloth and wipe all external surfaces of the device. If the surface has become contaminated with proteinaceous material, remove with a light detergent solution before surfacing cleaning with alcohol.

Prevent any fluid from entering the device.

WARNING 	Electric Shock hazard: Do not immerse the device.
WARNING 	Do NOT attempt to sterilise the device
WARNING 	Do NOT use a 100% alcohol or any other solvent to clean the device as this may cause damage to the casing surface and display.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

WARNING 	Protect the device from ingress of liquid. Should any liquid enter the device, discontinue use immediately and refer to an authorised service agent
---	---

8. SERVICE INFORMATION

8.1 SERVICE INDICATOR Service intervals are based on 1000 hours of pump operating time measured when the unit is providing suction. This equates to approximately 3000 aspiration cycles before the pump requires servicing.

The Service Indicator LED (arrowed) will remain illuminated indicating the pump now requires servicing. After 50 hours of usage without servicing, the Service Indicator will flash continuously indicating servicing is now overdue.

Illumination of the service indicator does not inhibit pump function.

Once illuminated, the pump must be returned for servicing as soon as possible.

WARNING 	No user serviceable parts inside.
---	---

8.2 RETURNING THE PUMP FOR SERVICE All devices to be returned must be prepared as described below for the protection of the servicing team and for safety during transport.

- Surface clean the pump as described in the Section above.
- Seal in a plastic bag and seal within a second plastic bag.
- **Place in the original packaging.**
- Enclose the following information:
 - Contact name
 - Centre address
 - Decontamination Certificate
 - Description of the fault or service required
 - Accompanying Order to authorise servicing – contact your local Customer Services Team for details.

WARNING 	No user serviceable parts inside.
	Do NOT include used consumables as these pose a significant contamination risk
	A decontamination certificated MUST be included with every returned pump. Repair or servicing cannot be commenced unless the service agent is in possession of this certificate.

8.3 SERVICE AGENTS

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

UK & Ireland:
IVF Synergy Ltd.
1st Floor Suite
Maritime House
Discovery Quay
Falmouth TR11 3XA
UNITED KINGDOM
Tel: +44 (0) 1872 487224
Email: service@ivfsynergy.com

Australian Service
Agent :
Rocket Medical Pty Ltd
Suite 6,
157 Gordon Street,
Port Macquarie NSW
2444
AUSTRALIA
Email:
aus@rocketmedical.com

Rest of World:
Charles Austen Pumps
Charles Austen Pumps Ltd
Royston Road
Byfleet, KT14 7NY
UNITED KINGDOM
Tel: +44 (0) 1932 355277
Email:
[Sales@charlesausten.com](mailto:sales@charlesausten.com)

UK Customer Services:
Rocket Medical plc.
Sedling Road
Washington NE38 9BZ
UNITED KINGDOM
Tel: +44 (0) 191 419
6988
Fax: +44 (0) 191 419
6989
Email:
customerservices@rocketmedical.com

9. WARRANTY

Oocyte Aspiration Pumps are sold by **Rocket Medical plc** under the warranties set forth in the following paragraphs. Such warranties are extended only with respect to the purchase of the Products directly from **Rocket Medical plc** as new merchandise and are extended to the first Buyer thereof, other than for resale.

For a period of TWENTY-FOUR (24) months from the date of shipment the Products are warranted to be free from functional defects in materials and workmanship and to conform to the description of the Products contained in the operating manual and accompanying labels, provided the same is properly operated under conditions of normal use, that annual maintenance and service is performed at an authorised **Rocket Medical plc** service facility

Removal of any QC seal voids the warranty.

The foregoing warranties shall not apply if the Products have been repaired other than by **Rocket Medical plc** or other than in accordance with written instructions provided by **Rocket Medical plc**, or altered by anyone other than **Rocket Medical plc**, or if the Products have been subject to misuse, negligence, or accident.

Rocket Medical plc 's sole and exclusive obligation and Buyer's sole and exclusive remedy under the above warranties is limited to repairing or replacing, free of charge, at **Rocket Medical plc** 's option, Products, which are reported to **Rocket Medical plc** by mail, telephone or email and which, if so advised by **Rocket® Medical plc**, is thereafter returned with a statement of the observed deficiency, not later than seven (7) days after the expiration date of the warranty, to **Rocket Medical plc** during normal business address, transport charges prepaid and which, upon **Rocket Medical plc**'s examination, is not found to conform with the above warranties.

Rocket Medical plc shall not be otherwise liable for any damages including but not limited to incidental damages, consequential damages or special damages.

There are no express or implied warranties which extend beyond the warranties herein above set forth. **Rocket Medical plc** makes no warranty of merchantability or fitness for a particular purpose with respect to the Products or parts thereof.

10. TECHNICAL SPECIFICATIONS

10.1 CLASSIFICATION

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

IEC 60601-1

Type of protection against electric shock: Class I.

Degree of protection against electric shock: Type B.

Vacuum type: high vacuum/low volume.

Suitable for continuous operation.

Not suitable for use in the presence of flammable gases.

Not suitable for use in conditions which expose the device to the ingress of water.


Not suitable for sterilisation.

10.2 SPECIFICATIONS

Power Input to Pump: 12VDC

Universal Power Supply: 100V to 240 VAC @ 50 / 60Hz Maximum current: 2.5A @ 12V

Service interval: 1000hrs of pump operation.:

	Dimensions: □ W - 248mm □ H - 86mm □ D - 194mm	Weight: • Unit – 2.56 Kg • Foot Switch:0.51 Kg
---	---	--

Suction Ranges:

User Set: -50mmHg to -300mmHg in -5mmHg increments

Pre-set: -500mmHg – nominal

Accurate to ±5% of full scale

10.3 EMC TABLES:

Table 1: Guidance and Manufacturer's Declaration – Electromagnetic Emissions

Emissions Test	Compliance	Electromagnetic environment - guidance
RF Emissions CISPR 11	Group 1	The Oocyte Aspiration Pump generates RF signals for its internal function. Therefore, its RF emissions are very low and not likely to cause any interference in nearby electronic equipment.
RF CISPR 11	Class A	The Oocyte Aspiration Pump is suitable for use in all establishments, other than domestic and those connected directly to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

Table 2: Guidance and Manufacturer's Declaration - Electromagnetic Immunity

The Oocyte Aspiration Pump (the device) is intended for use within the electromagnetic environment specified below. The operator of a device should ensure that it is used within such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic environment Guidance

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

Electrostatic discharge (ESD) IEC 61000-4-2	+/- 6kV contact +/- 8kV air	+/- 6kV contact +/- 8kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient burst IEC 61000-4-4	+/- 2kV for power supply lines	+/- 2kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	+/- 1kV differential mode +/- 2kV common mode	+/- 1kV differential mode +/- 2kV common mode	Mains power quality should be that of a typical commercial or hospital environment
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 610004-11	<5% Ut (>95% dip in Ut) for 0.5 cycle 40% Ut (60% dip in Ut) for 5 cycles 70% Ut (30% dip in Ut) for 25 cycles <5% UT (>95% dip in Ut) for 5 seconds	<5% Ut (>95% dip in Ut) for 0.5 cycles 40% Ut (60% dip in Ut) for 5 cycles 70% Ut (30% dip in Ut) for 25 cycles <5% UT (>95% dip in Ut) for 5 seconds	Mains power quality should be that a typical commercial or hospital environment. If the operator of the device requires continued operation during mains interruptions, it is recommended that the device be powered from an uninterruptable power supply.
Note: Ut is the a.c. mains voltage prior to application of the test level.			

Table 3: Guidance and Manufacturer's Declaration - Electromagnetic Immunity

The device is intended for use within the electromagnetic environment specified below. The customer or the operator of a should ensure that it is used within such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic environment Guidance
			Portable and mobile RF communications equipment should be no closer to any part of the PSU device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter
			Recommended Separation Distance
Conducted Rf	3 Vrms	3Vrms	$D = 1.2 \sqrt{P}$
IEC 61000-4-6	150 kHz to 80MHz		

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

Radiated RF	3V/m	3V/m	$D = 1.2 \sqrt{P}$ 80 MHz to 800 MHz
IEC 61000-4-3	80 MHz to 2.5 GHz		$D = 1.2 \sqrt{P}$ 800 MHz to 300 GHz
			Where P is the maximum power output [power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m) Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey a, should be less than the compliance level in each frequency in each frequency range.
			Interference may occur in the vicinity of equipment marked with the following symbol
Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.			
Note2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and /or people.			
a. Field strengths from fixed transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device. b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Table 4: Recommended Separation Distances Between Portable and Mobile RF Communications Equipment and the Oocyte Aspiration Pump.

The device is intended for use in the electromagnetic environment specified below. The operator of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.			
	Separation Distance According to Frequency of Transmitter m		
Read Maximum Output Power of Transmitter	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
W	$d = 1.2 \sqrt{P}$	$d = 1.2 \sqrt{P}$	$d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	3.8	3.8	7.3
10	3.8	3.8	7.3
100	12	12	23

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

For transmitters rated at a maximum output power not listed above, the recommended separation distances d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and/or people.

II. DISPOSAL

This device, its accessories and the consumables used with it, should be handled and disposed of in accordance with policy of the healthcare setting and with regard to all applicable regulations, including but without limitation to, those pertaining to human health & safety and care of the environment. Particular care to be taken with sharp objects such as needles due to the risk of infection, cuts and environmental hazards.

Waste Electrical and Electronic Equipment: In Europe, disposal of this device must be undertaken with regard to the WEEE directive (2002/96/EC). Do not dispose of this device in municipal waste.

In the event a serious incident related to this device occurs, the event should be reported to Rocket Medical at pncf@rocketmedical.com as well as to the competent health authority in the country that the user/patient resides.

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