

CRAFT™ Suction Pumps



CRAFT™ Suction and CRAFT™ DUO-VAC™ Suction Pumps have been the gold standard for reliable oocyte aspiration for over 30 years. They have been continuously developed to provide smooth, low-volume vacuum at pre-determined negative pressures, so allowing simple, low-turbulence, oocyte recovery.

- HIGHLY RELIABLE:** Rocket Medical CRAFT™ pumps have a proven 30-year history of reliability and excellent value.
- SIMPLE, SAFE OPERATION** with air-operated foot switches on a solid moulded base. Lightweight connection set across the range to improve ease of handling and reliability.
- LARGE, CLEAR GAUGE** with increased resolution to improve setting accuracy.
- ULTRA-QUIET** diaphragm pump makes all Rocket Medical pumps easy to live with in the operating room.
- FULLY VARIABLE VACUUM**
Medium Vacuum: 0 to -200mmHg. High Vacuum: -440mmHg.

DESCRIPTION	TYPE	ORDER CODE
CRAFT™ Suction Pump Complete with air-operated vacuum footswitch, 5 single patient use water trap sets, mains supply cord, User Manual.	220-240V AC 50Hz	R29655
	110-120V AC 60Hz	R29654
CRAFT™ DUO-VAC™ Suction Pump Complete with air-operated high/medium vacuum footswitch, 5 single patient use water trap sets, mains supply cord, User Manual.	220-240V AC 50Hz	R29660
	110-120V AC 60Hz	R29661
Water Trap Set for CRAFT™ Suction Pumps 2.0m patient connection tube with male luer connector. Supplied sterile, for single use in cartons of 10 units.		R57685

